

**ADPT001**

**Pyrogen free Test Tubes (10 X 75 mm)**

For the testing / detection of endotoxin content by LAL (Limulus Amebocyte Lysate) method

### **Composition**

Disposable Borosilicate Glass Test Tubes with Plain end are made of borosilicate glass to reduce pH changes and contaminants potentially leached from soda lime glass.

### **Appearance:**

Shrink wrap pack contains 5 packs of each 10 number of 10 x 75 mm flint glass tubes tightly wrapped in aluminium foil

### **Principle:**

Each pack of 10 number of 10 X 75 mm test tubes are tightly wrapped in aluminium foil and subjected to Depyrogenation for 1 hour at 250°C.

### **Endotoxin Content:**

< 0.005 EU / Tube

