

ADPT002

Pyrogen free Dilution Tubes (16 X 100 mm)

For the testing / detection of endotoxin content by LAL (Limulus Amebocyte Lysate) method

Composition

Disposable Borosilicate Glass Test Tubes with Plain end are made of borosilicate glass to reduce pH changes and contaminants potentially leached from soda lime glass.

Appearance:

Shrink wrap pack contains 4 packs of each 5 number of 16 x 100 mm flint glass tubes tightly wrapped in aluminium foil

Principle:

Each pack of 5 number of 16 X 100 mm test tubes are tightly wrapped in aluminium foil and subjected to Depyrogenation for 1 hour at 250°C.

Endotoxin Content:

< 0.005 EU / Tube

